



examination table

7,135 m

33° 58' 48" N, 76° 1' 18"



**EVEREST**  
veterinary technology

# KUN: EXAMINATION TABLE

## PRACTICE, ECONOMIC AND DETACHABLE



7,135 m  
33° 58' 48" N, 76° 1' 18"

### INNOVATIVE AND ECONOMICAL

*Examination table which, due to the quality of materials and easy mounting system, maintains the quality levels of the furniture designed and manufactured by EVEREST, with an economical price.*

### IT COME DISASSEMBLED

*The KUN table comes disassembled, in a flat pack. This helps us to offer you a high quality product at reduced costs and less risk of incidences of transport.*

### FEATURES

#### Materials:

100% Stainless steel AISI (18/10)  
+ wood under the surface to avoid metallic noises that scare the animal.

SURFACE: 112 X 55 cm

HEIGHT: 85 cm

CAPACITY: 200 kg

MEASURING LEGS: 50 x 50 mm

ADJUSTABLE LEG: Compensates for possible ground levels.

2 MEDICAL RAILS:

Stainless steel bars, Universal format for attachment of accessories.

**MR040105**

### PRACTICE AND DETACHABLE

Specially designed for veterinary use.  
Manufactured in stainless steel AISI 304 (18/10)  
Fast and safe assembly. Easy to transport.

### SOUNDPROOF

Anti-noise Surface: Wood base that absorbs metallic sounds scare the animal.

### ERGONOMIC DESIGN

Ergonomic surface, without bottom bars.  
Modern design with straight lines that facilitate cleaning.

### STABLE

Adjustable feet that compensate for possible ground level.

### 2 UNIVERSAL RAILS

In stainless steel to place any type of accessory.