

# Anaesthesia Equipment

## **McKinley** type 2 Anaesthetic Machines

- Steel trolley made with anticorrosive protection and epoxy painting.
- 75 mm high-resistance plastic castors with individual front brakes.
- Holds 1 x F size cylinder secured by belt and clamp, minimal lifting required.
- Flowmeter system: 0-4 LMP for O<sub>2</sub>.
- Protection bar for flowmeter valves
- Selectatec system for vapouriser mounting.
- Stainless steel absorber valves
- Quick and simple absorber dismounting system
- 1-year warranted
- Common front gas outlet standart for all circuits.
- CO<sub>2</sub> autoclavable for 0.4 Kg of soda lime. Built-in APL valve in the upper part of the absorber.
- Oxygen flush max 65 lpm (more than enough for animals up to 80Kg)



Optional: McKinley type 8 built-in volume ventilator electronic and digital.

1 SODA  
LIME

+

5 TRACHEAL  
TUBE

+

T D'AYRES  
+ 0,5 liter  
rebreathing bag

+

BAIN  
+ 2 liter  
rebreathing bag

+

CIRCULAR CIRCUIT  
+ 1,2 and 3 liter  
rebreathing bag

OPTIONAL



VAPORISER ISOFLURANE

VV010...



O<sub>2</sub> CONCENTRATOR

CO050003

[www.everest-tecnovet.com](http://www.everest-tecnovet.com)

[info@everest-tecnovet.com](mailto:info@everest-tecnovet.com) +34 936 680 6715

 **EVEREST**  
tecnologia veterinaria

# Anaesthesia Equipment

**McKinley**  
Anaesthetic Machines ONE

## ULTRA COMPACT AND ECONOMICAL



· FULLY AUTOCLAVABLE  
PLASTIC CANISTER



· CLOSED STAINLESS STEEL  
SUPPORT THAT ENSURES  
SECURITY IN THE INTERNAL  
CONNECTIONS



· VAPORISER INTEGRATED,  
WITHOUT CONNECTIONS.  
MAXIMUM SECURITY




· EXCLUSIVE EQUIPMENT WITH  
INTEGRATED STRIP 5 SOCKETS



Using the latest tools of industrial design, Everest has created an ultra-compact anaesthesia equipment that includes all the elements that are part of larger equipments, reducing space to a minimum, increasing the comfort of use.

Thus, we present the new McKinley One, with a lot of possibilities and an affordable price, to adapt to the needs of all veterinary clinics, without losing the maximum security that requires an anaesthesia machine.

1 SODALIME + 5 TRACHEAL + T D'AYRES + CIRCULAR CIRCUIT  
TUBE + 0,5 liter rebreathing bag + 1,2 and 3 liter rebreathing bag



MCKINLEY ONE  
MOBILE STAND

EQUIPMENT  
+ RODABLE SUPPORT  
+ CONCENTRATOR  
+ ISOFLURANE VAPORISER



MCKINLEY ONE  
TABLE TOP / WALL MOUNTED

EQUIPMENT  
+ ISOFLURANE VAPORISER

\*Options: Table top, wall mounted or mobile stand, the new McKinley ONE will suit any surgery.  
\*Implementation by tutorial.

[www.everest-tecnovet.com](http://www.everest-tecnovet.com)

[info@everest-tecnovet.com](mailto:info@everest-tecnovet.com) +34 936 680 6715

**EVEREST**  
tecnologia veterinaria