



examination table

7,135 m

33° 58' 48" N, 76° 1' 18"



EVEREST
tecnovet

You Tube

NUN: EXAMINATION TABLE FOR THE MOST DEMANDING VETERINARIANS



UNBEATABLE

NUN
examination table

7,135 m
33° 58' 48" N, 76° 1' 18"

FEATURES:

Materials:
100% Stainless steel AISI (18/10) + wood under the surface to avoid metallic noises that scare the animal.

SURFACE: 120 x 60 cm

HEIGHT: 85 cm

CAPACITY: 200 kg

MEASURING LEGS: 40 x 40 mm and 30 x 30 mm reinforcements of stainless pipe.

ADJUSTABLE LEG: Compensates for possible ground levels.

2 MEDICAL RAILS:
Stainless steel bars, Universal format for attachment of accessories.

One of the oldest tables manufactured by Everest. It actually exists since before our founding. Until today, many of our customers still have their tables purchased in the 1990s. And they enjoy them.

RESISTANT AND DURABLE

Specially designed for veterinary use. Manufactured in stainless steel AISI 304 (18/10). Single piece, without screws or assemblies.

SOUNDPROOF

Anti-noise Surface: Wood base that absorbs metallic sounds scare the animal.

ERGONOMIC DESIGN

Central Bottom Bar: Avoid hitting your feet while you work.

STABLE

Adjustable feet that compensate for possible ground level.

2 UNIVERSAL RAILS

In stainless steel to place any type of accessory.



MR040103