



 **McKinley** type2  
Anaesthetic Machines

▶ 6.194 m

📍 63° 4' N 151° 0' O



# EQUIPO DE ANESTESIA: McKinley type 2



- Soporte portátil de acero inoxidable con estantes para monitores.
- Ruedas de 75 mm de alta resistencia con frenos individuales.
- Soporte rodable del equipo y de la botella de O<sub>2</sub> (10 l) de acero pintado.
- Circuito circular de metal cromado.
- Caudalímetro de precisión de oxígeno 0-4 LPM para O<sub>2</sub>.
- Barra de protección para la válvula del caudalímetro.
- Anclaje de vaporizador tipo Selectatec.
- Circuito circular con absorbedor de CO<sub>2</sub> de fácil uso con canister de 400 gr. (con 400 g de cal sodada que ayuda a mantener muy bajos los flujos de oxígeno en cualquier tamaño de paciente).
- Válvulas de acero inoxidable.
- Rápido y simple sistema de cambio para la cal sodada.

**McKinley** type 2  
Anaesthetic Machines

▶ 6.194 m

📍 63° 4' N 151° 0' O

## OPCIONALES:

VAPORIZADOR ISOFLURANO  
VV010...



CONCENTRADOR O<sub>2</sub>  
C0050001



## INCLUYE:

- 1 CAL SODADA
- + 5 TUBOS TRAQUEALES
- + T D'AYRES + BALÓN DE 0,5 L
- + CIRCUITO CERRADO
- + 3 BALONES DE 1, 2, 3 LITROS

**MCKTWO**

CONSULTE OFERTAS VIGENTES



 **McKinley** type2  
Anaesthetic Machines

▶ 6.194 m

📍 63° 4' N 151° 0' O



# ANESTHESIA MACHINE: McKinley type 2



- Steel trolley made with anticorrosive protection and epoxy painting.
- 75 mm high-resistance plastic castors with individual front brakes.
- Holds 1 x F size cylinder secured by belt and clamp, minimal lifting required.
- Flowmeter system: 0-4 LMP for O<sub>2</sub>.
- Protection bar for flowmeter valves.
- Selectatec system for vaporiser mounting.
- Stainless steel absorber valves.
- Common front gas outlet standart for all circuits.
- CO<sub>2</sub> autoclavable for 0.4 Kg of soda lime. Built-in APL valve in the upper part of the absorber.
- Oxygen flush max 65 lpm (more than enough for animals up to 80Kg).

**McKinley** type 2  
Anaesthetic Machines

▶ 6.194 m

📍 63° 4' N 151° 0' O

## OPTIONAL:

ISOFLURANE VAPORIZER  
VV010...



OXIGEN CONCENTRATOR  
C0050003



## INCLUDE:

- + 1 SODA LIME
- + 5 TRACHEAL TUBES
- + T D'AYRES
- + 0,5 liter rebreathing bag
- + BAIN
- + 2 liter rebreathing bag
- + CIRCULAR CIRCUIT
- + 1, 2, 3 liter rebreathing bag

**MCKTWO**

ASK FOR CURRENT OFFERS