



mesa de examen

7,135 m

33° 58' 48" N, 76° 1' 18"



EVEREST  
tecnovet

You Tube

# NUN: LA MESA DE EXAMEN PARA LOS VETERINARIOS MÁS EXIGENTES

*Una de las mesas más antiguas fabricada por Everest. En realidad existe desde antes de nuestra fundación. Hasta hoy, muchos de nuestros clientes aún tienen sus mesas compradas en los 90. Y las disfrutan.*



## INMEJORABLE

### RESISTENTE Y DURABLE

Diseñada especialmente para uso veterinario.  
Fabricada en acero inoxidable AISI 304 (18/10).  
Pieza única, sin tornillos ni montajes.

### INSONORIZADA

Superficie anti-ruídos: Base de madera que absorbe los sonidos metálicos que asusten al animal.

### DISEÑO ERGONÓMICO

Barra inferior central: Evita golpearse los pies mientras trabaja.

### ESTABLE

Patas ajustables que compensan posibles desniveles del suelo.

### 2 RAÍLES UNIVERSALES

En acero inoxidable para colocar cualquier tipo de accesorio.



# NUN

## mesa de examen

7,135 m

33° 58' 48" N, 76° 1' 18"

### CARACTERÍSTICAS

#### MATERIALES:

100% Acero inoxidable AISI (18/10)  
+ Madera bajo la superficie para evitar ruidos metálicos que asusten al animal.

SUPERFICIE: 120 x 60 cm

ALTURA: 85 cm

ALTURA REGULABLE: 80-11 cm

CAPACIDAD: 200 Kg

Medida patas: 40 x 40 mm y refuerzos de 30 x 30 mm de tubo inoxidable.

PATA AJUSTABLE: Compensa posibles desniveles del suelo.

#### 2 RAÍLES MÉDICOS:

Barras en acero inox, formato universal para sujeción de accesorios.



**MR040101**

ALTURA FIJA  
85 cm

**MR040104**

ALTURA REGULABLE  
80-110 cm

[www.everest-tecnovet.com](http://www.everest-tecnovet.com)

+34 93 668 6715 · [info@everest-tecnovet.com](mailto:info@everest-tecnovet.com)

 **EVEREST**  
tecnología veterinaria



examination table

▶ 7,135 m

📍 33° 58' 48" N, 76° 1' 18"



EVEREST  
tecnovet

You Tube

# NUN: EXAMINATION TABLE

## FOR THE MOST DEMANDING VETERINARIANS

*One of the oldest tables manufactured by Everest. It actually exists since before our founding. Until today, many of our customers still have their tables purchased in the 1990s. And they enjoy them.*



### RESISTANT AND DURABLE

Specially designed for veterinary use. Manufactured in stainless steel AISI 304 (18/10). Single piece, without screws or assemblies.

### SOUNDPROOF

Anti-noise Surface: Wood base that absorbs metallic sounds scare the animal.

### ERGONOMIC DESIGN

Central Bottom Bar: Avoid hitting your feet while you work.

### STABLE

Adjustable feet that compensate for possible ground level.

### 2 UNIVERSAL RAILS

In stainless steel to place any type of accessory.



# NUN

examination table

7,135 m  
33° 58' 48" N, 76° 1' 18"

### FEATURES:

#### Materials:

100% Stainless steel AISI (18/10) + wood under the surface to avoid metallic noises that scare the animal.

SURFACE: 120 x 60 cm

FIXED HEIGHT: 85 cm

ADJUSTABLE HEIGHT: 80-110 cm

CAPACITY: 200 kg

MEASURING LEGS: 40 x 40 mm and 30 x 30 mm reinforcements of stainless pipe.

ADJUSTABLE LEG: Compensates for possible ground levels.

#### 2 MEDICAL RAILS:

Stainless steel bars, Universal format for attachment of accessories.



### MR040101

FIXED HEIGHT  
85 cm

### MR040104

ADJUSTABLE HEIGHT  
80-110 cm