

CIRCUITOS RESPIRATORIOS



ÚNICO CIRCUITO EXCLUSIVO PARA USO VETERINARIO

EVEREST ofrece una gama exclusiva de circuitos anestésicos para uso exclusivo veterinario, disponibles en diferentes configuraciones según peso del animal, para anestesia, quirófano y sala de recuperación.

CONEXIONES FÁCILES Y SEGURAS

Conectores hembra de 22 mm, fabricados en material suave, ligero y flexible, garantizan un ajuste seguro y facilitan la conexión/desconexión al equipo de anestesia.


FABRICACIÓN LIMPIA

Todos los circuitos se fabrican en salas blancas.

CARACTERÍSTICAS

- Diferentes tamaños: para animales de <10Kg, 10 - 30 Kg y >30 Kg
- Tubos de 120 cm de PVC de ánima lisa.
- Conectores en forma de "Y" y de codo
- CON y SIN trampa de agua.

Todos los productos tienen la marca CE y cumplen con las regulaciones ISO 13485.



"Y" Ø10 mm	"Y" Ø15 mm	"Y" Ø19 mm	"Y" Ø15 mm trampa de agua	"Y" Ø19 mm trampa de agua
Conector "Soft" Ø10/22 mm Tubo de Ánima lisa Ø10 mm Longitud 120 cm Peso paciente <10 Kg	Conector "Soft" Ø15/22 mm Tubo de Ánima lisa Ø15 mm Longitud 120 cm Peso paciente 10-30 Kg	Conector "Soft" Ø 22 mm Tubo de Ánima lisa Ø19 mm Longitud 120 cm Peso paciente >30 Kg	Conector "Soft" Ø15/22 mm Tubo de Ánima lisa Ø15 mm Longitud 120 cm Peso paciente 10-30 Kg	Conector "Soft" Ø22 mm Tubo de Ánima lisa Ø19 mm Longitud 120 cm Peso paciente ≥30 Kg
CA101001	CA101501	CA102201	CA101505	CA102205

BREATHING CIRCUITS



THE ONLY CIRCUIT EXCLUSIVELY FOR VETERINARY USE

EVEREST offers an exclusive range of anesthetic circuits for veterinary use only, available in different combinations according to animal weight, for anesthesia, surgery and recovery room.

EASY AND SAFE CONNECTIONS

22 mm female connectors made of soft, light and flexible material, which guarantee a secure fit and facilitate connection/disconnection to the anesthesia machine.


CLEAN MANUFACTURING

All circuits are manufactured in clean rooms.

FEATURES

- Different sizes: for animals <10Kg, 10 - 30 Kg and >30 Kg.
- 120 cm smooth-bore PVC tubing.
- "Y" and elbow connectors.
- WITH and WITHOUT water trap.

All products are CE marked and comply with ISO 13485 regulations.



"Y" Ø10 mm	"Y" Ø15 mm	"Y" Ø19 mm	"Y" Ø15 mm with water trap	"Y" Ø19 mm with water trap
"Soft" Connector Ø10/22 mm Smoothbore tube Ø10 mm Length 120 cm Patient weight <10 Kg	"Soft" Connector Ø15/22 mm Smoothbore tube Ø15 mm Length 120 cm Patient weight 10-30 Kg	"Soft" Connector Ø 22 mm Smoothbore tube Ø19 mm Length 120 cm Patient weight >30 Kg	"Soft" Connector Ø15/22 mm Smoothbore tube Ø15 mm Length 120 cm Patient weight 10-30 Kg	"Soft" Connector Ø22 mm Smoothbore tube Ø19 mm Length 120 cm Patient weight ≥30 Kg
CA101001	CA101501	CA102201	CA101505	CA102205